


# Spectroscopy Application Note

## Analysis of Cobalt Base Alloys

Cobalt base alloys are typically associated with applications requiring wear, corrosion, or heat resistance. The chemical composition of these alloys signifies the basis for which application they are best suited.

LECO manufactures three glow discharge optical emission spectrometers designed for optimum elemental determination in ferrous and nonferrous materials. A "glow discharge" source uniformly removes material from the sample surface. The separation of sampling from sample excitation reduces the effect of metallurgical and chemical history inherent in all samples. The excitation of primarily ground state atom lines gives rise to less complex spectra, minimizing and/or eliminating many interferences. Calibration curves are linear and cover a very wide dynamic range. The same wavelengths are often used for both high and low concentration ranges.



### Typical Sample Results

<b>Stellite 1 #LO 2228</b>											
	<b>Cr</b>	<b>W</b>	<b>Mo</b>	<b>C</b>	<b>Fe</b>	<b>Ni</b>	<b>Si</b>	<b>Mn</b>	<b>B</b>	<b>Cu</b>	<b>Nb</b>
Avg (%)	25.73	11.11	0.36	2.22	1.64	1.42	1.09	0.11	0.0074	0.03	0.0086
Std Dev	0.078	0.12	0.0022	0.022	0.016	0.014	0.0044	0.0018	0.0001	0.0002	0.0004
RSD (%)	0.3	1.12	0.6	1.0	0.96	1.0	0.4	1.59	0.97	0.59	5.12
<b>Stellite 21 #LO 2234</b>											
	<b>Cr</b>	<b>W</b>	<b>Mo</b>	<b>C</b>	<b>Fe</b>	<b>Ni</b>	<b>Si</b>	<b>Mn</b>	<b>B</b>	<b>Cu</b>	<b>Nb</b>
Avg (%)	26.81	0.41	5.5	0.22	0.69	2.48	0.63	0.52	N/A	0.063	N/A
Std Dev	0.098	0.0042	0.0017	0.0002	0.0057	0.008	0.0005	0.0069	—	0.0065	—
RSD (%)	0.37	1.03	0.03	0.11	0.83	0.32	0.08	1.32	—	10	—
<b>Alloy 800 #LO 2240</b>											
	<b>Cr</b>	<b>W</b>	<b>Mo</b>	<b>C</b>	<b>Fe</b>	<b>Ni</b>	<b>Si</b>	<b>Mn</b>	<b>B</b>	<b>Cu</b>	<b>Nb</b>
Avg (%)	17.55	0.021	28.49	0.01	1.34	1.33	3.11	0.03	0.008	N/A	N/A
Std Dev	0.11	0.0031	0.067	0.00004	0.0059	0.016	0.012	0.0003	0.00001	—	—
RSD (%)	0.62	15.0	0.23	0.37	0.44	1.23	0.37	0.88	0.14	—	—
<b>CO-Cr-W #LO 2238</b>											
	<b>Cr</b>	<b>W</b>	<b>Mo</b>	<b>C</b>	<b>Fe</b>	<b>Ni</b>	<b>Si</b>	<b>Mn</b>	<b>B</b>	<b>Cu</b>	<b>Nb</b>
Avg (%)	29.32	19.66	0.25	0.63	2.77	5.44	0.79	0.26	0.023	0.024	0.061
Std Dev	0.025	0.066	0.0013	0.0049	0.017	0.01	0.0019	0.0013	0.0002	0.00005	0.0001
RSD (%)	0.09	0.34	0.52	0.77	0.61	0.18	0.24	0.48	0.74	0.21	0.19
<b>UNCO 50 #LO 2520</b>											
	<b>Cr</b>	<b>W</b>	<b>Mo</b>	<b>C</b>	<b>Fe</b>	<b>Ni</b>	<b>Si</b>	<b>Mn</b>	<b>B</b>	<b>Cu</b>	<b>Nb</b>
Avg (%)	27.54	1.06	0.011	0.14	28.3	0.61	0.77	0.84	0.0006	0.016	0.20
Std Dev	0.059	0.022	0.0002	0.0004	0.082	0.0005	0.0024	0.0053	0.00001	0.00004	0.0012
RSD (%)	0.21	2.0	1.48	0.3	0.29	0.074	0.31	0.63	1.67	0.26	6.18

**GDS-Series**

## Typical Results—Certified Cobalt Standards

Brammer BS 170B, Cobalt Alloy 6B											
	Cr	W	Mo	C	Fe	Ni	Si	Mn	B	Cu	Nb
Avg (%)	30.51	3.95	0.81	1.12	0.92	2.50	0.71	1.40	0.00041	0.016	0.027
Std Dev	0.011	0.011	0.0034	0.023	0.0058	0.012	0.0039	0.0024	0.00002	0.0002	0.00083
RSD (%)	0.037	0.28	0.42	2.1	0.63	0.49	0.54	0.17	6	1.3	3.1
Certified	30.52	3.98	0.81	1.12	0.91	2.49	0.71	1.41	0.00040	0.016	0.027

IARM 96A Stellite 25											
	Cr	W	Mo	C	Fe	Ni	Si	Mn	B	Cu	Nb
Avg (%)	20.11	15.16	0.28	0.010	2.29	10.53	0.17	1.52	0.0030	0.040	0.069
Std Dev	0.0010	0.023	0.00064	0.00038	0.0050	0.017	0.0012	0.0045	0.00001	0.00011	0.00086
RSD (%)	0.0052	0.15	0.23	3.8	0.22	0.16	0.69	0.3	0.56	0.300	1.2
Certified	20.17	15.20	0.28	0.010	2.29	10.53	0.17	1.52	0.0030	0.040	0.070



  

MBH 14942-L Stellite 8											
	Cr	W	Mo	C	Fe	Ni	Si	Mn	B	Cu	Nb
Avg (%)	29.01	1.58	6.28	0.097	1.07	0.27	1.09	0.48	N/A	N/A	N/A
Std Dev	0.023	0.0046	0.0098	0.0002	0.0026	0.0008	0.0020	0.0009	—	—	—
RSD (%)	0.079	0.29	0.16	0.17	0.25	0.30	0.18	0.19	—	—	—
Certified	28.56	1.60	6.39	0.10	1.07	0.27	1.07	0.49	—	—	—

Analyses shown were run on the LECO GDS-400A.

## Calibration

Linear working curves are found with LECO glow discharge spectrometers. Single wavelength lines often cover full concentration ranges. Linear calibrations correlate to low reference material consumption and few spectral interferences.



## Sample Preparation

The techniques specified in ASTM Practices and Standards apply to glow discharge. As a non-thermal sputtering source, glow discharge does not rely on high temperatures to melt and volatilize samples. Craftsman-like machining techniques are avoided, and surface preparation with belt or disc grinders will easily provide a uniform surface sufficient for glow discharge. Cobalt base alloys are typically abraded with a 120 grit zirconium oxide belt/disc.

## Accessories

Sample surface preparation (LECO BG-30, VP-160 or other suitable equipment).

## Calibration Standards

Based on customer requirements; NIST, Brammer, MBH, ARMI or other suitable standards.

## Typical Analysis Times

Start-up and Pre-burn	~30 seconds
Analyze	~15 seconds
Total	~45 seconds



LECO Corporation • 3000 Lakeview Ave. • St. Joseph, MI 49085-2396  
 Phone: 800-292-6141 • Fax: 269-982-8977  
 info@leco.com • www.leco.com • ISO-9001 No. FM 24045

LECO is a registered trademark of LECO Corporation  
 6/03-REV1