

# Inorganic Application Note



## Carbon and Sulfur in Lime

### **Instrument**

CS-200, CS-300, CS-400 and CS-444 Series Carbon and Sulfur Determinators

### **Calibration Standard**

NIST or other suitable standards  
Suggestions: LECO 502-319 @ 1.36% C and 1.53% S,  
NIST SRM 2690 Fly Ash @ 0.15% S

### **Accessories**

LECO 528-018 Ceramic Crucibles (preheated),  
LECO 773-579 Metal Scoop, 763-266 LECOCEL or  
763-263 LECOCEL III, and 502-231 or 501-077 Iron Chip  
Accelerators (for carbon-only instruments substitute 501-263 Copper Accelerator for LECOCEL or LECOCEL III)

### **Sample Weight**

~0.20 to 0.25 grams

### **Sample Preparation**

Sample should be uniform powder.

### **Instrument Settings**

Power Level:	maximum
Pre-analyze Purge*:	5 seconds
Pre-analyze Delay*:	10 seconds
Carbon Minimum Time-out:	50 seconds
Carbon Comparator Level:	1.00%
Sulfur Minimum Time-out:	60 seconds
Sulfur Comparator Level:	1.00%

*\*For improved precision at low carbon levels (<0.1%), a pre-analyze purge of 10 and delay of 25 seconds is recommended.*

### **Note**

Samples containing water of hydration (crystallization) or hydroxides could experience diminished sulfur recovery.



# Carbon/Sulfur

## Method

1. Preheat ceramic crucibles in a muffle or tube furnace at 1200 to 1300°C for not less than 15 minutes or at 950 to 1050°C for not less than 40 minutes. Remove the crucibles from the furnace, cool for 1 to 2 minutes, and place in a desiccator for storage. If the crucibles are not used within four hours, they should be rebaked.
2. Determine the blank.
  - a. Enter 1.000 g weight into weight stack.
  - b. Add one level 773-579 Metal Scoop of LECOCEL followed by one level metal scoop of iron chip accelerator to a preheated crucible.
  - c. Place crucible on furnace pedestal and analyze.
  - d. Repeat steps 2a through 2c a minimum of five times.
  - e. Enter blank following routine outlined in operator's instruction manual.
3. Calibrate.
  - a. Weigh ~0.20 to 0.25 g calibration standard into the center of the preheated crucible. Enter the weight into the weight stack.
  - b. Add one level 773-579 Metal Scoop of LECOCEL followed by one level metal scoop of iron chip accelerator covering the calibration sample.
  - c. Place crucible on furnace pedestal and analyze.
  - d. Repeat steps 3a through 3c a minimum of five times and calibrate the instrument following the calibration procedure as outlined in the operator's instruction manual.
  - e. Verify the calibration by analyzing the calibration standard again. It should fall within the expected tolerance. If not repeat steps 3a through 3e.
4. Analyze samples.
  - a. Weigh ~0.20 to 0.25 g sample into the center of the crucible. Enter the weight into the weight stack.
  - b. Add one level 773-579 Metal Scoop of LECOCEL followed by one level metal scoop of iron chip accelerator covering the sample.
  - c. Place crucible on furnace pedestal and analyze.

## Typical Results

	Weight (g)	Carbon (%)	Sulfur (%)		Weight (g)	Carbon (%)	Sulfur (%)
LECO 502-319	0.2056	1.36	1.53	Quicklime	0.1951	1.914	0.0168
Ore	0.2035	1.36	1.53	Sample A	0.2120	1.842	0.0174
@1.36% C,	0.2024	1.36	1.53		0.2120	1.910	0.0161
1.53% S	0.2175	1.36	1.54				
	0.2130	1.36	1.52	Quicklime	0.2020	0.606	0.0220
				Sample B	0.1993	0.653	0.0213
NIST 2690	0.1966	0.393	0.151		0.2072	0.595	0.0218
Fly Ash	0.2072	0.392	0.150				
@ 0.15% S	0.2000	0.386	0.150	Dolomitic Lime	0.2171	0.209	0.0249
Carbon is	0.2019	0.383	0.149	Sample A	0.2067	0.235	0.0240
not certified	0.1960	0.383	0.151		0.1951	0.213	0.0232
				Dolomitic Lime	0.2217	0.0952	0.0149
				Sample B	0.2040	0.0918	0.0151
					0.2173	0.0914	0.0152



**LECO Corporation**  
 3000 Lakeview Ave. • St. Joseph, MI 49085-2396  
 Phone: 800-292-6141 • Fax: 269-982-8977  
 info@leco.com • www.leco.com • ISO-9001 No. FM 24045