

Organic Application Note



Sulfur in Oils

Accessories

529-204 Crucibles, 502-301 Tungsten Trioxide or 502-321 COM-CAT accelerator, eye droppers or disposable pipettes

Sample Weight

0.09 to 0.12 g

Calibration Standard

NIST or other LECO suitable material

Furnace Temperature

1350°C

Analysis Time

60 to 300 seconds (Sulfur)

Comparator Level

1% (Sulfur)

Notes

- If the sample soots (indicated by a darkening of the anhydron in the primary anhydron tube), then a high bias of the results is possible. To overcome sooting, add ~3 g of 501-426 COM-AID on top of the ~1 g of Tungsten Trioxide to cover the sample.
- 502-321 COM-CAT can be substituted for 502-301 Tungsten Trioxide.

Procedure

1. Prepare instrument as outlined in the operator manual.
2. Condition the instrument by analyzing a few coal samples, until the readings are stable.
3. Calibrate the system by analyzing at least three replicates of the standard used, using the following procedure.
 - a. Using an eyedropper, weigh ~0.1 g of sample into a 529-204 Crucible containing ~1 g 502-301 Tungsten Trioxide or 502-321 COM-CAT.
 - b. Click the mouse on Analyze. When the "load furnace" message appears on screen, slide the sample into the combustion tube.
4. Analyze unknown samples following steps 3a and 3b.



S-144DR

Typical Results

Sample	Weight (g)	% Sulfur
501-099 LECO	0.1137	4.03
Paraffin Oil	0.1159	4.04
	0.1020	4.02
	0.1049	4.07
	0.1165	4.03
	0.1185	4.10
	0.1138	4.00
	0.1137	3.98
	0.1167	4.04
	0.1189	4.06
	Mean	= 4.04
	Std. Dev.	= 0.03



LECO Corporation

3000 Lakeview Ave. • St. Joseph, MI 49085-2396

Phone: 800-292-6141 • Fax: 269-982-8977

info@leco.com • www.leco.com • ISO-9001 No. FM 24045