

Organic Application Note

Nitrogen in Oil/Amines

Accessories 502-186 Tin Foil Cups, 501-614 Spatula, Disposable Eyedroppers, 501-439 Paraffin Oil, 501-426 COM-AID

Sample Weight 0.1 g

Calibration Standard 502-092 EDTA, Glycine Solution (see reverse side for detailed instructions), or other suitable standard

Furnace Temperature 950°C

Flow Profile All High

Atmospheric Blank (N) 0.04 (verified by analyzing paraffin oil with COM-AID)

Crucible Changing Interval ~100 analyses using 614-961-110 Porous Crucible

Analysis Time ~170 seconds

Procedure

1. Prepare the instrument by following the procedure as outlined in the operator's instruction manual (i.e. check gas supplies, perform any required maintenance, perform leak checks, etc.).
2. Analyze blanks (gas) until a plateau is reached. Analyze three to five additional blanks and set blank using these data.
3. Analyze five EDTA standards at 0.2 g and drift correct (if using the PC option). *NOTE: Each method on PC requires prior calibration with multiple weights of EDTA (0.035 to 0.4 g). If PC is not installed, analyze five EDTA standards and calibrate using the DSP screen menu.*
4. Weigh ~0.1 g oil into a 502-186 Tin Foil Cup containing ~0.1 g 601-427 COM-AID, add 0.4 to 0.5 g additional COM-AID, seal tin foil to avoid trapping air or losing sample, and analyze.
5. Analyze a standard at end of set to verify calibration.

Typical Results

Sample	Weight (g)	% Nitrogen
Amine #1	0.0982	1.00
	0.1074	0.98
	0.1074	1.00
	Average =	0.99
	Std. Dev. =	0.014
Amine #2	0.0901	0.143
	0.0951	0.131
	0.0876	0.141
	Average =	0.138
	Std. Dev. =	0.007

Sample	Weight (g)	% Nitrogen
Amine #3	0.1072	4.80
	0.0931	4.82
	0.1030	4.84
	0.1054	4.83
	0.1079	4.79
	0.1107	4.79
	Average =	4.81
	Std. Dev. =	0.021
Lube Oil	0.1041	0.101
	0.1059	0.111
	0.1050	0.100
	0.1054	0.100
	Average =	0.103
Std. Dev. =	0.005	



FP-528

GLYCINE SOLUTION PREPARATION

1. The following formula can be used to make a specific concentration:

$$G = \frac{1}{[0.18658 - 0.01] C}$$

where: C = desired nitrogen concentration as percent
G = grams of glycine powder

Example for 1% solution:

$$G = \frac{1}{[0.18658 - 0.01] 1.00} = 5.663$$

NOTE: A quick reference chart, shown below, shows the grams of glycine powder needed to reach given concentrations.

- Place a flask on the balance and tare. The flask should be large enough to hold 100 ml (where 100 g = 100 ml).
- Add the amount of glycine calculated in step 1 and record the weight.
- Tare the balance and add 100 g of distilled water and record the weight.
- Seal the flask and mix the contents.
- To figure the exact concentration:

$$\% \text{ Nitrogen} = \frac{G (18.658)}{G + W}$$

where: G = weight in grams of glycine recorded in step 3
W = weight in grams of water recorded in step 4

- If the distilled water is not pure, determining the nitrogen concentration may be necessary.
 - Using the liquid injector, analyze five samples of distilled water.
 - Average the nitrogen content of the five samples (A).
 - Add this average to % nitrogen calculated for the calibration solution.

example: To make a calibration solution of approximately 0.3% nitrogen:

$$\begin{aligned} \text{where: } G &= 1.672 \text{ g} \\ W &= 99.824 \text{ g} \\ A &= 0.004\% \end{aligned}$$

$$\frac{1.672(18.654)}{(1.672 + 99.824)} + 0.004 = 0.311\%$$

QUICK REFERENCE CONCENTRATION TABLE

Nitrogen Concentration Grams of Glycine

0.10%	0.539 g
0.30%	1.634 g
0.50%	2.754 g
0.75%	4.188 g
1.00%	5.663 g



LECO Corporation • 3000 Lakeview Ave. • St. Joseph, MI 49085-2396
Phone: 800-292-6141 • Fax: 269-982-8977
info@leco.com • www.leco.com • ISO-9001 No. FM 24045