

# Organic Application Note

## Total Mercury in Industrial Wastes

### Accessories

614-822-102 Small Nickel Boats.

Note: Boats should be pre-baked at 400° C or analyzed (without a sample) before loading a sample.

**Sample Weight** 50 to 250 mg (0.050 to 0.250 g)

### Calibration Standard

NIST SRM 1633b (coal fly ash), NIST SRM 2709 (soil), and NIST SRM 2711 (soil); or equivalent Certified Reference Materials.

**Furnace Temperature** 550° C for catalyst; 750° C for decomp

**Analysis Time** ~ 6 minutes

### Method Profile

Drying Time: 60 seconds

Decomposition Time: 250 seconds

Cuvette Cleaning Cycle: 45 seconds

Peak Used: Low for wastewater; high for others

NOTE: Method for Quicksilver Windows® Software Version 2.0.

### Procedure

- Determine the blank as follows:
  - Enter "Blank" from the drop-down menu under the "Name" column.
  - Click "Analyze", the door will open and the nickel loop will be presented.
  - Carefully place a 614-822-102 Small Nickel Boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window, the door will close and the analysis sequence will start automatically.
  - Repeat steps 1a through 1d two more times. The system and boats will be purged of any interfering elements.
- Calibrate the instrument as defined in the instructional manual:
  - Analyze various sample weights of a relevant reference material in accordance to the absolute amount of mercury required to calibrate an appropriate dynamic range. The calibration samples are weighed into the 614-822-102 Small Nickel Boat.
  - Enter each calibration sample with the appropriate ID code from the drop-down menu, and sample weight from an external balance measurement.
  - Click "Analyze", the door will open and the nickel loop will be presented.
  - If there is a boat in the nickel loop, remove it and keep for later use.
  - Carefully place the calibration sample boat into the nickel loop using clean tweezers.
  - Click "OK" in the "Load Sample" window, the door will close and the analysis sequence will start automatically.
  - Repeat steps 2a through 2f as per the calibration procedures.  
Note: The first analyzed sample after a long delay should be discarded. This sample should be considered a conditioner for the system, and not used for the actual calibration.
  - Complete a calibration by following the calibration procedure as outlined in the manual.
  - Verify the calibration by analyzing one of the calibration samples again. It should be within the expected tolerances. If not, repeat steps 2a through 2i.



**AMA254**

3. Analyze the samples as follows:
  - a. Weigh ~60 mg of the high concentration sample into a 614-822-102 Small Nickel Boat.  
*NOTE: Use ~200 mg for low concentration samples.*
  - b. Enter a sample identification in the Name column and the sample weight in the Mass column.
  - c. Click "Analyze", the door will open and the nickel loop will be presented.
  - d. If there is a boat in the nickel loop, remove it and keep for later use.
  - e. Carefully place the sample boat into the nickel loop using clean tweezers.
  - f. Click "OK" in the "Load Sample" window, the door will close and the analysis sequence will start automatically.

### Typical Results

Sample	Weight(g)	Hg (ppm)	Sample	Weight(g)	Hg (ppm)
Oil/Sludge	0.0498	1.61	Wastewater	0.1563	0.004
	0.0619	1.69		0.2235	0.003
	0.0558	1.69		0.2042	0.003
	<b>Avg. (ppm)</b>	<b>1.66</b>		<b>Avg. (ppm)</b>	<b>0.0033</b>
Sand	0.0604	4.95	Shale	0.1781	0.025
	0.0566	5.08		0.1711	0.022
	0.0538	4.98		0.1671	0.022
	0.0534	5.08		<b>Avg. (ppm)</b>	<b>0.023</b>
<b>Avg. (ppm)</b>	<b>5.02</b>				



LECO Corporation • 3000 Lakeview Ave. • St. Joseph, MI 49085-2396  
 Phone: 800-292-6141 • Fax: 269-982-8977  
 info@leco.com • www.leco.com • ISO-9001 No. FM 24045